

## Nuclear Matrix Protein p84 Rabbit mAb

Catalog No: #59524

Package Size: #59524-1 50ul #59524-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Nuclear Matrix Protein p84 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	Nuclear Matrix Protein p84 Antibody detects endogenous levels of total Nuclear Matrix Protein p84
Immunogen Description	A synthesized peptide derived from human Nuclear Matrix Protein p84
Other Names	Death domain containing protein p84N5; HPR1; hTREX84; Nuclear matrix protein p84; p84N5; THO complex 1; Tho1; Thoc1;
Accession No.	Uniprot:Q96FV9
Uniprot	Q96FV9
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200

## Images



Western blot analysis of Nuclear Matrix Protein p84 expression in HepG2 cell lysate.

## Product Description

Regulates transcriptional elongation of a subset of genes. Participates in an apoptotic pathway which is characterized by activation of caspase-6, increases in the expression of BAK1 and BCL2L1 and activation of NF-kappa-B.

## Background

---

Regulates transcriptional elongation of a subset of genes. Participates in an apoptotic pathway which is characterized by activation of caspase-6, increases in the expression of BAK1 and BCL2L1 and activation of NF-kappa-B.

---

Note: This product is for in vitro research use only