

## Follistatin Rabbit mAb

Catalog No: #59538

Package Size: #59538-1 50ul #59538-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Follistatin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Rat
Specificity	Follistatin Antibody detects endogenous levels of total Follistatin
Immunogen Description	A synthesized peptide derived from human Follistatin
Other Names	Activin binding protein; Follistatin; FS; fst; Human follistatin gene;
Accession No.	Uniprot:P19883
Uniprot	P19883
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000

## Images



Western blot analysis of Follistatin expression in HepG2 cell lysate.

## Product Description

Binds directly to activin and functions as an activin antagonist. Specific inhibitor of the biosynthesis and secretion of pituitary follicle stimulating hormone (FSH).

## Background

---

Binds directly to activin and functions as an activin antagonist. Specific inhibitor of the biosynthesis and secretion of pituitary follicle stimulating hormone (FSH).

---

Note: This product is for in vitro research use only