

Cyclophilin F Rabbit mAb

Catalog No: #59539

Package Size: #59539-1 50ul #59539-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Cyclophilin F Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	Cyclophilin F Antibody detects endogenous levels of total Cyclophilin F
Immunogen Description	A synthesized peptide derived from human Cyclophilin F
Other Names	Cyclophilin 3; cyclophilin D; Cyclophilin F; Cyp D; CyP M; CYP3; CypD; hCyP3; Mitochondrial cyclophilin; PPlase; PPlase F; Ppif; Rotamase;
Accession No.	Uniprot:P30405
Uniprot	P30405
Calculated MW	18kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50

Images



Western blot analysis of Cyclophilin F expression in HeLa cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Product Description

PPIases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.

Background

PPIases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.

Note: This product is for in vitro research use only