

GABARAPL2 Rabbit mAb

Catalog No: #59552

Package Size: #59552-1 50ul #59552-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GABARAPL2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	GABARAPL2 Antibody detects endogenous levels of total GABARAPL2
Immunogen Description	A synthesized peptide derived from human GABARAPL2
Other Names	ATG8C; FLC3A; GABA(A) receptor-associated protein-like 2; Gabarapl2; Ganglioside expression factor 2; GATE16; GEF2; General protein transport factor p16;
Accession No.	Uniprot:P60520
Uniprot	P60520
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of GABARAPL2 expression in (1) HeLa cell lysate; (2) Mouse spleen lysate; (3) Rat brain lysate.

Product Description

Involved in intra-Golgi traffic. Modulates intra-Golgi transport through coupling between NSF activity and SNAREs activation. It first stimulates the ATPase activity of NSF which in turn stimulates the association with GOSR1.

Background

Involved in intra-Golgi traffic. Modulates intra-Golgi transport through coupling between NSF activity and SNAREs activation. It first stimulates the ATPase activity of NSF which in turn stimulates the association with GOSR1.

Note: This product is for in vitro research use only