

PMP22 Rabbit mAb

Catalog No: #59556

Package Size: #59556-1 50ul #59556-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PMP22 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Rat
Specificity	PMP22 Antibody detects endogenous levels of total PMP22
Immunogen Description	A synthesized peptide derived from human PMP22
Other Names	CMT1A; CMT1E; DSS; GAS3; Growth Arrest Specific 3; HMSNIA; HNPP; Peripheral myelin protein 22; PMP22; Sp110; Trembler;
Accession No.	Uniprot:Q01453
Uniprot	Q01453
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of PMP22 expression in 293T cell lysate.

Product Description

Involved in growth regulation, and in myelination in the peripheral nervous system.

Background

Involved in growth regulation, and in myelination in the peripheral nervous system.

Note: This product is for in vitro research use only