

Calretinin Rabbit mAb

Catalog No: #59564

Package Size: #59564-1 50ul #59564-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Calretinin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC FC
Species Reactivity	Human Mouse Rat
Specificity	Calretinin Antibody detects endogenous levels of total Calretinin
Immunogen Description	A synthesized peptide derived from human Calretinin
Other Names	29 kDa calbindin; CAB29; CAL 2; CAL2; CALB 2; CALB2; Calbindin 2 29kDa; Calbindin 2; Calbindin D29K; Calbindin2; Calretinin; CR;
Accession No.	Uniprot:P22676
Uniprot	P22676
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:1000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of Calretinin expression in SH-SY5Y cell lysate.

Product Description

Calretinin is a calcium-binding protein which is abundant in auditory neurons.

Background

Calretinin is a calcium-binding protein which is abundant in auditory neurons.

Note: This product is for in vitro research use only