

NDRG1 Rabbit mAb

Catalog No: #59565

Package Size: #59565-1 50ul #59565-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NDRG1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	NDRG1 Antibody detects endogenous levels of total NDRG1
Immunogen Description	A synthesized peptide derived from human NDRG1
Other Names	Protein NDRG1; Differentiation-related gene 1 protein; DRG-1; RTP; Rit42; NDRG1; CAP43; DRG1; RTP; targ1; TDD5; tdds;
Accession No.	Uniprot:Q92597
Uniprot	Q92597
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of NDRG1 expression in HeLa cell lysate.

Product Description

May have a growth inhibitory role.

Background

May have a growth inhibitory role.

Note: This product is for in vitro research use only