

StAR Rabbit mAb

Catalog No: #59567

Package Size: #59567-1 50ul #59567-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	StAR Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	StAR Antibody detects endogenous levels of total StAR
Immunogen Description	A synthesized peptide derived from human StAR
Other Names	StAR; StARD1; START domain containing 1; Steroid acute regulatory protein;
Accession No.	Uniprot:P49675
Uniprot	P49675
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:50~1:200

Images



Western blot analysis of StAR expression in SK-BR-3 cell lysate.

Product Description

Steroidogenic acute regulatory protein (StAR) plays a significant role in cholesterol transport from the cytoplasmic outer membrane to the inner mitochondrial membrane.

Background

Steroidogenic acute regulatory protein (StAR) plays a significant role in cholesterol transport from the cytoplasmic outer membrane to the inner mitochondrial membrane.

Note: This product is for in vitro research use only