

ETS1 Rabbit mAb

Catalog No: #59571

Package Size: #59571-1 50ul #59571-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ETS1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse
Specificity	ETS1 Antibody detects endogenous levels of total ETS1
Immunogen Description	A synthesized peptide derived from human ETS1
Other Names	C ets 1 protein; ETS 1; ETS1; ETS1 oncogene; EWSR2; Oncogene ETS1; P54; Protein C-ets-1;
Accession No.	Uniprot:P14921
Uniprot	P14921
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 FC 1:80

Images



Western blot analysis of Insulin Receptor R expression in Daudi cell lysate.

Product Description

Transcription factor.

Background

Transcription factor.

Note: This product is for in vitro research use only