

## Thymidylate Synthase Rabbit mAb

Catalog No: #59572

Package Size: #59572-1 50ul #59572-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Thymidylate Synthase Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC IP FC
Species Reactivity	Human Mouse Rat
Specificity	Thymidylate Synthase Antibody detects endogenous levels of total Thymidylate Synthase
Immunogen Description	A synthesized peptide derived from human Thymidylate Synthase
Other Names	d TMP synthase; HsT422; Thymidylate synthase; TMS; TS; TSase; Tyms; TYMS protein; Tyms thymidylate synthetase;
Accession No.	Uniprot:P04818
Uniprot	P04818
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

## Images



Western blot analysis of Thymidylate Synthase expression in HeLa cell lysate.

## Product Description

Contributes to the de novo mitochondrial thymidylate biosynthesis pathway.

## Background

---

Contributes to the de novo mitochondrial thymidylate biosynthesis pathway.

---

Note: This product is for in vitro research use only