

PKC beta 1 Rabbit mAb

Catalog No: #59578

Package Size: #59578-1 50ul #59578-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PKC beta 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	PKC beta 1 Antibody detects endogenous levels of total PKC beta 1
Immunogen Description	A synthesized peptide derived from human PKC beta 1
Other Names	PKC beta; PKCB; Prkcb; PRKCB I; PRKCB1; PRKCB2; Protein kinase C beta;
Accession No.	Uniprot:P05771
Uniprot	P05771
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:100

Images



Western blot analysis of PKC beta 1 expression in (1) HeLa cell lysate; (2) NIH/3T3 cell lysate.



Immunohistochemical analysis of paraffin-embedded mouse colon, using PKC beta 1 Antibody.

Product Description

Plays a key role in B-cell activation and function by regulating BCR-induced NF-kappa-B activation and B-cell survival.

Background

Plays a key role in B-cell activation and function by regulating BCR-induced NF-kappa-B activation and B-cell survival.

Note: This product is for in vitro research use only