

USP22 Rabbit mAb

Catalog No: #59588

Package Size: #59588-1 50ul #59588-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	USP22 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human
Specificity	USP22 Antibody detects endogenous levels of total USP22
Immunogen Description	A synthesized peptide derived from human USP22
Other Names	Deubiquitinating enzyme 22; Ubiquitin carboxyl terminal hydrolase 22; Usp22; USP3L;
Accession No.	Uniprot:Q9UPT9
Uniprot	Q9UPT9
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

Images



Western blot analysis of USP22 expression in HeLa cell lysate.

Product Description

Histone deubiquitinating component of the transcription regulatory histone acetylation (HAT) complex SAGA. Catalyzes the deubiquitination of both histones H2A and H2B, thereby acting as a coactivator.

Background

Histone deubiquitinating component of the transcription regulatory histone acetylation (HAT) complex SAGA. Catalyzes the deubiquitination of both histones H2A and H2B, thereby acting as a coactivator.

Note: This product is for in vitro research use only