

## Apg12 Rabbit mAb

Catalog No: #59601

Package Size: #59601-1 50ul #59601-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Apg12 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	Apg12 Antibody detects endogenous levels of total Apg12
Immunogen Description	A synthesized peptide derived from human Apg12
Other Names	APG12-like; APG12L; ATG12; Autophagy 12; FBR93; HAPG12;
Accession No.	Uniprot:O94817
Uniprot	O94817
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200

## Images



Western blot analysis of Apg12 expression in HepG2 cell lysate.

## Product Description

Ubiquitin-like protein required for autophagy. Conjugated to ATG3 and ATG5.

## Background

Ubiquitin-like protein required for autophagy. Conjugated to ATG3 and ATG5.

---

Note: This product is for in vitro research use only