

## Huntingtin Rabbit mAb

Catalog No: #59607

Package Size: #59607-1 50ul #59607-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Huntingtin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	Huntingtin Antibody detects endogenous levels of total Huntingtin
Immunogen Description	A synthesized peptide derived from human Huntingtin
Other Names	HD; HDH; HTT; Huntingtin; IT15; ZHD;
Accession No.	Uniprot:P42858
Uniprot	P42858
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:100~1:500 FC 1:100

## Images



Western blot analysis of Huntingtin expression in HeLa cell lysate.

## Product Description

May play a role in microtubule-mediated transport or vesicle function.

## Background

May play a role in microtubule-mediated transport or vesicle function.

---

Note: This product is for in vitro research use only