

RANTES Rabbit mAb

Catalog No: #59622

Package Size: #59622-1 50ul #59622-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	RANTES Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human
Specificity	RANTES Antibody detects endogenous levels of total RANTES
Immunogen Description	A synthesized peptide derived from human RANTES
Other Names	SISd; eoCP; SCYA5; RANTES; TCP228; D17S136E; SIS-delta;
Accession No.	Uniprot:P13501
Uniprot	P13501
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 ICC/IF 1:50~1:200

Images



Western blot analysis of RANTES expression in RANTES recombinant protein lysate.

Product Description

Chemoattractant for blood monocytes, memory T-helper cells and eosinophils. Causes the release of histamine from basophils and activates eosinophils. Binds to CCR1, CCR3, CCR4 and CCR5.

Background

Chemoattractant for blood monocytes, memory T-helper cells and eosinophils. Causes the release of histamine from basophils and activates eosinophils. Binds to CCR1, CCR3, CCR4 and CCR5.

Note: This product is for in vitro research use only