

MYBBP1A Rabbit mAb

Catalog No: #59627

Package Size: #59627-1 50ul #59627-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MYBBP1A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	MYBBP1A Antibody detects endogenous levels of total MYBBP1A
Immunogen Description	A synthesized peptide derived from human MYBBP1A
Other Names	cb486; Mybbp1a; P160; p160MBP; p67MBP; PAP2; RP23 48A2.3;
Accession No.	Uniprot:Q9BQG0
Uniprot	Q9BQG0
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:1000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of MYBBP1A expression in HEK293 cell lysate.

Product Description

May activate or repress transcription via interactions with sequence specific DNA-binding proteins. Repression may be mediated at least in part by histone deacetylase activity (HDAC activity).

Background

May activate or repress transcription via interactions with sequence specific DNA-binding proteins. Repression may be mediated at least in part by histone deacetylase activity (HDAC activity).

Note: This product is for in vitro research use only