

GLP2 Rabbit mAb

Catalog No: #59634

Package Size: #59634-1 50ul #59634-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GLP2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	GLP2 Antibody detects endogenous levels of total GLP2
Immunogen Description	A synthesized peptide derived from human GLP2
Other Names	GCG; GLP 2; GRPP;
Accession No.	Uniprot:P01275
Uniprot	P01275
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of GLP2 expression in GLP2 recombinant protein lysate.

Product Description

Human GLP2 (glucagon like peptide 2) is a peptide hormone of the glucagon family produced by the L cells of the intestinal mucosa and originates from the same prohormone as glucagon.

Background

Human GLP2 (glucagon like peptide 2) is a peptide hormone of the glucagon family produced by the L cells of the intestinal mucosa and originates from the same prohormone as glucagon.

Note: This product is for in vitro research use only