

## GCLM Rabbit mAb

Catalog No: #59636

Package Size: #59636-1 50ul #59636-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	GCLM Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	GCLM Antibody detects endogenous levels of total GCLM
Immunogen Description	A synthesized peptide derived from human GCLM
Other Names	GCLM; GCS light chain; GLCLR; Glutamate cysteine ligase regulatory subunit;
Accession No.	Uniprot:P48507
Uniprot	P48507
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:100

## Images



Western blot analysis of GCLM expression in K562 cell lysate.

## Product Description

Sulfur metabolism; glutathione biosynthesis; glutathione from L-cysteine and L-glutamate: step 1/2.

## Background

Sulfur metabolism; glutathione biosynthesis; glutathione from L-cysteine and L-glutamate: step 1/2.

---

Note: This product is for in vitro research use only