

## PGC1 beta Rabbit mAb

Catalog No: #59637

Package Size: #59637-1 50ul #59637-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PGC1 beta Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	PGC1 beta Antibody detects endogenous levels of total PGC1 beta
Immunogen Description	A synthesized peptide derived from human PGC1 beta
Other Names	PERC; PGC1; PPARGC1; Ppargc1b;
Accession No.	Uniprot:Q86YN6
Uniprot	Q86YN6
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IP 1:50

## Images



Western blot analysis of PGC1 beta expression in (1) HeLa cell lysate; (2) Mouse spleen lysate; (3) Rat spleen lysate.

## Product Description

Plays a role of stimulator of transcription factors and nuclear receptors activities. Activates transcriptional activity of estrogen receptor alpha, nuclear respiratory factor 1 (NRF1) and glucocorticoid receptor in the presence of glucocorticoids.

## Background

---

Plays a role of stimulator of transcription factors and nuclear receptors activities. Activates transcriptional activity of estrogen receptor alpha, nuclear respiratory factor 1 (NRF1) and glucocorticoid receptor in the presence of glucocorticoids.

---

Note: This product is for in vitro research use only