

## Aldolase Rabbit mAb

Catalog No: #59640

Package Size: #59640-1 50ul #59640-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Aldolase Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	Aldolase Antibody detects endogenous levels of total Aldolase
Immunogen Description	A synthesized peptide derived from human Aldolase
Other Names	ALDA; Aldo1; ALDOA; Aldolase 1; Aldolase A; GSD12; HEL S 87p;
Accession No.	Uniprot:P04075
Uniprot	P04075
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

## Images



Western blot analysis of Aldolase expression in A549 cell lysate.

## Product Description

Plays a key role in glycolysis and gluconeogenesis. In addition, may also function as scaffolding protein.

## Background

Plays a key role in glycolysis and gluconeogenesis. In addition, may also function as scaffolding protein.

---

Note: This product is for in vitro research use only