

TPOR Rabbit mAb

Catalog No: #59653

Package Size: #59653-1 50ul #59653-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	TPOR Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	TPOR Antibody detects endogenous levels of total TPOR
Immunogen Description	A synthesized peptide derived from human TPOR
Other Names	C MPL; CD110; mpl; MPLV; THCYT2; Thrombopoietin receptor; TPO R; TPOR;
Accession No.	Uniprot:P40238
Uniprot	P40238
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of TPOR expression in (1)A431 cell lysate; (2)NIH/3T3 cell lysate.

Product Description

Receptor for thrombopoietin. May represent a regulatory molecule specific for TPO-R-dependent immune responses.

Background

Receptor for thrombopoietin. May represent a regulatory molecule specific for TPO-R-dependent immune responses.

Note: This product is for in vitro research use only