

Oct-2 Rabbit mAb

Catalog No: #59654

Package Size: #59654-1 50ul #59654-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Oct-2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	Oct-2 Antibody detects endogenous levels of total Oct-2
Immunogen Description	A synthesized peptide derived from human Oct-2
Other Names	NF A2; Oct 2; Oct-2; OTF-2; OTF2; POU2F2;
Accession No.	Uniprot:P09086
Uniprot	P09086
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:50 FC 1:100

Images



Western blot analysis of Oct-2 expression in Ramos cell lysate.

Product Description

Transcription factor that specifically binds to the octamer motif (5'-ATTTGCAT-3'). Regulates transcription in a number of tissues in addition to activating immunoglobulin gene expression.

Background

Transcription factor that specifically binds to the octamer motif (5'-ATTTGCAT-3'). Regulates transcription in a number of tissues in addition to activating immunoglobulin gene expression.

Note: This product is for in vitro research use only