

## MMP17 Rabbit mAb

Catalog No: #59660

Package Size: #59660-1 50ul #59660-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	MMP17 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	MMP17 Antibody detects endogenous levels of total MMP17
Immunogen Description	A synthesized peptide derived from human MMP17
Other Names	MMP 17; Mmp17; MT MMP4; MT-MMP 4; MT4 MMP; MT4-MMP; MT4MMP; MTMMP4;
Accession No.	Uniprot:Q9ULZ9
Uniprot	Q9ULZ9
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:1000 IHC 1:50~1:200 ICC/IF 1:50~1:200

## Images



Western blot analysis of MMP17 expression in HeLa cell lysate.

## Product Description

Endopeptidase that degrades various components of the extracellular matrix, such as fibrin. May be involved in the activation of membrane-bound precursors of growth factors or inflammatory mediators, such as tumor necrosis factor-alpha. May also be involved in tumoral process.

## Background

---

Endopeptidase that degrades various components of the extracellular matrix, such as fibrin. May be involved in the activation of membrane-bound precursors of growth factors or inflammatory mediators, such as tumor necrosis factor-alpha. May also be involved in tumoral process.

---

Note: This product is for in vitro research use only