

NR1D1 Rabbit mAb

Catalog No: #59672

Package Size: #59672-1 50ul #59672-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NR1D1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	NR1D1 Antibody detects endogenous levels of total NR1D1
Immunogen Description	A synthesized peptide derived from human NR1D1
Other Names	EAR1; hRev; Nr1d1; Rev erbAalpha; Reverba; THRA1; THRAL;
Accession No.	Uniprot:P20393
Uniprot	P20393
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:20

Images



Western blot analysis of NR1D1 expression in HeLa cell lysate.

Product Description

Functions as a constitutive transcriptional repressor. Possible receptor for triiodothyronine.

Background

Functions as a constitutive transcriptional repressor. Possible receptor for triiodothyronine.

Note: This product is for in vitro research use only