

Matrin 3 Rabbit mAb

Catalog No: #59677

Package Size: #59677-1 50ul #59677-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Matrin 3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	Matrin 3 Antibody detects endogenous levels of total Matrin 3
Immunogen Description	A synthesized peptide derived from human Matrin 3
Other Names	KIAA0723; Matr3; Matrin3; MPD2; VCPDM;
Accession No.	Uniprot:P43243
Uniprot	P43243
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:100~1:500 ICC/IF 1:50~1:200 FC 1:150

Images



Western blot analysis of Matrin 3 expression in (1) HEK293 cell lysate; (2) PC12 cell lysate.

Product Description

May play a role in transcription or may interact with other nuclear matrix proteins to form the internal fibrogranular network. In association with the SFPQ-NONO heteromer may play a role in nuclear retention of defective RNAs.

Background

May play a role in transcription or may interact with other nuclear matrix proteins to form the internal fibrogranular network. In association with the SFPQ-NONO heteromer may play a role in nuclear retention of defective RNAs.

Note: This product is for in vitro research use only