

ILF2 Rabbit mAb

Catalog No: #59678

Package Size: #59678-1 50ul #59678-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ILF2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP FC
Species Reactivity	Human Mouse Rat
Specificity	ILF2 Antibody detects endogenous levels of total ILF2
Immunogen Description	A synthesized peptide derived from human ILF2
Other Names	ILF 2; MGC8391; NF45;
Accession No.	Uniprot:Q12905
Uniprot	Q12905
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50 FC 1:100

Images



Western blot analysis of ILF2 expression in A431 cell lysate.

Product Description

Appears to function predominantly as a heterodimeric complex with ILF3. This complex may regulate transcription of the IL2 gene during T-cell activation. It can also promote the formation of stable DNA-dependent protein kinase holoenzyme complexes on DNA.

Background

Appears to function predominantly as a heterodimeric complex with ILF3. This complex may regulate transcription of the IL2 gene during T-cell activation. It can also promote the formation of stable DNA-dependent protein kinase holoenzyme complexes on DNA.

Note: This product is for in vitro research use only