

## EB3 Rabbit mAb

Catalog No: #59682

Package Size: #59682-1 50ul #59682-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	EB3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	EB3 Antibody detects endogenous levels of total EB3
Immunogen Description	A synthesized peptide derived from human EB3
Other Names	EB3; EBF3; MAPRE3; RP/EB Family; RP3;
Accession No.	Uniprot:Q9UPY8
Uniprot	Q9UPY8
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

## Images



Western blot analysis of EB3 expression in Ratmuscle cell lysate.

## Product Description

May be involved in microtubule polymerization, and spindle function by stabilizing microtubules and anchoring them at centrosomes. May play a role in cell migration.

## Background

---

May be involved in microtubule polymerization, and spindle function by stabilizing microtubules and anchoring them at centrosomes. May play a role in cell migration.

---

Note: This product is for in vitro research use only