

Uteroglobin Rabbit mAb

Catalog No: #59689

Package Size: #59689-1 50ul #59689-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Uteroglobin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	Uteroglobin Antibody detects endogenous levels of total Uteroglobin
Immunogen Description	A synthesized peptide derived from human Uteroglobin
Other Names	Blastokinin; CC10; CC16; CCPBP; CCSP; SCGB1A; UGB; UP1; Uteroglobin;
Accession No.	Uniprot:P11684
Uniprot	P11684
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 FC 1:50

Images



Western blot analysis of Uteroglobin expression in A549 cell lysate.

Product Description

Binds phosphatidylcholine, phosphatidylinositol, polychlorinated biphenyls (PCB) and weakly progesterone, potent inhibitor of phospholipase A2.

Background

Binds phosphatidylcholine, phosphatidylinositol, polychlorinated biphenyls (PCB) and weakly progesterone, potent inhibitor of phospholipase A2.

Note: This product is for in vitro research use only