

FNDC5 Rabbit mAb

Catalog No: #59693

Package Size: #59693-1 50ul #59693-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FNDC5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	FNDC5 Antibody detects endogenous levels of total FNDC5
Immunogen Description	A synthesized peptide derived from human FNDC5
Other Names	FNDC 5; FRCP2; irisin;
Accession No.	Uniprot:Q8NAU1
Uniprot	Q8NAU1
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of FNDC5 expression in Rat muscle cell lysate.

Product Description

Irisin: Contrary to mouse, may not be involved in the beneficial effects of muscular exercise, nor in the induction of browning of human white adipose tissue.

Background

Irisin: Contrary to mouse, may not be involved in the beneficial effects of muscular exercise, nor in the induction of browning of human white adipose tissue.

Note: This product is for in vitro research use only