

## Versican Rabbit mAb

Catalog No: #59696

Package Size: #59696-1 50ul #59696-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Versican Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	Versican Antibody detects endogenous levels of total Versican
Immunogen Description	A synthesized peptide derived from human Versican
Other Names	CSPG2; ERVR; GHAP; PGM; VCAN; Versican; WGN1;
Accession No.	Uniprot:P13611
Uniprot	P13611
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200

## Images



Western blot analysis of Versican expression in Mouse heart cell lysate.

## Product Description

May play a role in intercellular signaling and in connecting cells with the extracellular matrix. May take part in the regulation of cell motility, growth and differentiation. Binds hyaluronic acid.

## Background

---

May play a role in intercellular signaling and in connecting cells with the extracellular matrix. May take part in the regulation of cell motility, growth and differentiation. Binds hyaluronic acid.

---

Note: This product is for in vitro research use only