

hnRNP U Rabbit mAb

Catalog No: #59697

Package Size: #59697-1 50ul #59697-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	hnRNP U Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	hnRNP U Antibody detects endogenous levels of total hnRNP U
Immunogen Description	A synthesized peptide derived from human hnRNP U
Other Names	HNRNPU; HNRPU; p120; pp120; SAFA;
Accession No.	Uniprot:Q00839
Uniprot	Q00839
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:20

Images



Western blot analysis of hnRNP U expression in K562 cell lysate.

Product Description

Component of the CRD-mediated complex that promotes MYC mRNA stabilization. Binds to pre-mRNA. Has high affinity for scaffold-attached region (SAR) DNA. Binds to double- and single-stranded DNA and RNA.

Background

Component of the CRD-mediated complex that promotes MYC mRNA stabilization. Binds to pre-mRNA. Has high affinity for scaffold-attached region (SAR) DNA. Binds to double- and single-stranded DNA and RNA.

Note: This product is for in vitro research use only