

LRG1 Rabbit mAb

Catalog No: #59698

Package Size: #59698-1 50ul #59698-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	LRG1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	LRG1 Antibody detects endogenous levels of total LRG1
Immunogen Description	A synthesized peptide derived from human LRG1
Other Names	A2GL; HMFT1766; LRG; LRG1;
Accession No.	Uniprot:P02750
Uniprot	P02750
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50

Images



Western blot analysis of LRG1 expression in Human plasma lysate.

Product Description

Transforming growth factor beta receptor binding.

Background

Transforming growth factor beta receptor binding.

Note: This product is for in vitro research use only