

ALDH1L1 Rabbit mAb

Catalog No: #59706

Package Size: #59706-1 50ul #59706-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ALDH1L1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	ALDH1L1 Antibody detects endogenous levels of total ALDH1L1
Immunogen Description	A synthesized peptide derived from human ALDH1L1
Other Names	10 FTHFDH; aldh1l1; FDH; FTHFD;
Accession No.	Uniprot:O75891
Uniprot	O75891
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of ALDH1L1 expression in 293T cell lysate.

Product Description

Highly expressed in liver, pancreas and kidney.

Background

Highly expressed in liver, pancreas and kidney.

Note: This product is for in vitro research use only