

## PARP2 Rabbit mAb

Catalog No: #59707

Package Size: #59707-1 50ul #59707-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PARP2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	PARP2 Antibody detects endogenous levels of total PARP2
Immunogen Description	A synthesized peptide derived from human PARP2
Other Names	pADPRT-2; hPARP-2; ADPRT2; ADPRTL2;
Accession No.	Uniprot:Q9UGN5
Uniprot	Q9UGN5
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:30

## Images



Western blot analysis of PARP2 expression in Raji cell lysate.

## Product Description

Involved in the base excision repair (BER) pathway, by catalyzing the poly(ADP-ribosyl)ation of a limited number of acceptor proteins involved in chromatin architecture and in DNA metabolism.

## Background

---

Involved in the base excision repair (BER) pathway, by catalyzing the poly(ADP-ribosyl)ation of a limited number of acceptor proteins involved in chromatin architecture and in DNA metabolism.

---

Note: This product is for in vitro research use only