

ELMO1 Rabbit mAb

Catalog No: #59708

Package Size: #59708-1 50ul #59708-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ELMO1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	ELMO1 Antibody detects endogenous levels of total ELMO1
Immunogen Description	A synthesized peptide derived from human ELMO1
Other Names	CED12; Ced12 homolog 1; ELMO 1;
Accession No.	Uniprot:Q92556
Uniprot	Q92556
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

Images



Western blot analysis of ELMO1 expression in Jurkat cell lysate.

Product Description

Involved in cytoskeletal rearrangements required for phagocytosis of apoptotic cells and cell motility. Acts in association with DOCK1 and CRK.

Background

Involved in cytoskeletal rearrangements required for phagocytosis of apoptotic cells and cell motility. Acts in association with DOCK1 and CRK.

Note: This product is for in vitro research use only