

TRP1 Rabbit mAb

Catalog No: #59712

Package Size: #59712-1 50ul #59712-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TRP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	TRP1 Antibody detects endogenous levels of total TRP1
Immunogen Description	A synthesized peptide derived from human TRP1
Other Names	CAS2; CATB; GP75; OCA3; TRP1;TYRP; TYRP1; TYRRP;
Accession No.	Uniprot:P17643
Uniprot	P17643
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50

Images



Western blot analysis of TRP1 expression in Human melanoma lysate.

Product Description

Oxidation of 5,6-dihydroxyindole-2-carboxylic acid (DHICA) into indole-5,6-quinone-2-carboxylic acid. May regulate or influence the type of melanin synthesized.

Background

Oxidation of 5,6-dihydroxyindole-2-carboxylic acid (DHICA) into indole-5,6-quinone-2-carboxylic acid. May regulate or influence the type of melanin synthesized.

Note: This product is for in vitro research use only