

## RNA Helicase A Rabbit mAb

Catalog No: #59719

Package Size: #59719-1 50ul #59719-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	RNA Helicase A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	RNA Helicase A Antibody detects endogenous levels of total RNA Helicase A
Immunogen Description	A synthesized peptide derived from human RNA Helicase A
Other Names	DDX 9; DDX9; dhx9; Leukophysin; LKP; NDH2; NDHII; RHA;
Accession No.	Uniprot:Q08211
Uniprot	Q08211
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:20

## Images



Western blot analysis of RNA Helicase A expression in Jurkat cell lysate.

## Product Description

Unwinds double-stranded DNA and RNA in a 3' to 5' direction. Alteration of secondary structure may subsequently influence interactions with proteins or other nucleic acids. Functions as a transcriptional activator. Component of the CRD-mediated complex that promotes MYC mRNA stability.

## Background

---

Unwinds double-stranded DNA and RNA in a 3' to 5' direction. Alteration of secondary structure may subsequently influence interactions with proteins or other nucleic acids. Functions as a transcriptional activator. Component of the CRD-mediated complex that promotes MYC mRNA stability.

---

Note: This product is for in vitro research use only