

p40 Rabbit mAb

Catalog No: #59721

Package Size: #59721-1 50ul #59721-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	p40 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human
Specificity	p40 Antibody detects endogenous levels of total p40
Immunogen Description	A synthesized peptide derived from human p40
Other Names	p40; RAB9p40; RABEPK; RGD1310612; RP11 65N13.1; RP23-446N16.2;
Accession No.	Uniprot:Q7Z6M1
Uniprot	Q7Z6M1
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 FC 1:50

Images



Western blot analysis of p40 expression in HepG2 cell lysate.

Product Description

Rab9 effector required for endosome to trans-Golgi network (TGN) transport.

Background

Rab9 effector required for endosome to trans-Golgi network (TGN) transport.

Note: This product is for in vitro research use only