

Ikaros Rabbit mAb

Catalog No: #59723

Package Size: #59723-1 50ul #59723-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Ikaros Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC FC
Species Reactivity	Human Mouse Rat
Specificity	Ikaros Antibody detects endogenous levels of total Ikaros
Immunogen Description	A synthesized peptide derived from human Ikaros
Other Names	hIk 1; IK1; IKAROS; Ikzf1; LYF1; ZNFN1A1;
Accession No.	Uniprot:Q13422
Uniprot	Q13422
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:100~1:500 FC 1:60

Images



Western blot analysis of Ikaros expression in MOLT4 cell lysate.

Product Description

Transcription regulator of hematopoietic cell differentiation (PubMed:17934067). Binds gamma-satellite DNA (PubMed:17135265, PubMed:19141594). Plays a role in the development of lymphocytes, B- and T-cells.

Background

Transcription regulator of hematopoietic cell differentiation (PubMed:17934067). Binds gamma-satellite DNA (PubMed:17135265, PubMed:19141594). Plays a role in the development of lymphocytes, B- and T-cells.

Note: This product is for in vitro research use only