

P70 S6 Kinase beta Rabbit mAb

Catalog No: #59741

Package Size: #59741-1 50ul #59741-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	P70 S6 Kinase beta Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	P70 S6 Kinase beta Antibody detects endogenous levels of total P70 S6 Kinase beta
Immunogen Description	A synthesized peptide derived from human P70 S6 Kinase beta
Other Names	p70 beta; P70S6K2; p70S6Kb; Rps6kb2; S6K beta 2; S6K beta; S6K2; SRK; STK14B;
Accession No.	Uniprot:Q9UBS0
Uniprot	Q9UBS0
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of P70 S6 Kinase beta expression in K562 cell lysate.

Product Description

Phosphorylates specifically ribosomal protein S6.

Background

Phosphorylates specifically ribosomal protein S6.

Note: This product is for in vitro research use only