

KAT8 / MYST1 / MOF Rabbit mAb

Catalog No: #59754

Package Size: #59754-1 50ul #59754-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	KAT8 / MYST1 / MOF Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	KAT8 / MYST1 / MOF Antibody detects endogenous levels of total KAT8 / MYST1 / MOF
Immunogen Description	A synthesized peptide derived from human KAT8 / MYST1 / MOF
Other Names	hMOF; MOF; MOZ; MYST1; SAS2; SAS3; YBF2; YBF2/SAS3; ZC2HC8;
Accession No.	Uniprot:Q9H7Z6
Uniprot	Q9H7Z6
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:1000 IHC 1:50~1:200 ICC/IF 1:100~1:500 IP 1:100 FC 1:100

Images



Western blot analysis of KAT8 / MYST1 / MOF expression in HeLa cell lysate.

Product Description

Histone acetyltransferase which may be involved in transcriptional activation. May influence the function of ATM.

Background

Histone acetyltransferase which may be involved in transcriptional activation. May influence the function of ATM.

Note: This product is for in vitro research use only