

SULT2A1 Rabbit mAb

Catalog No: #59757

Package Size: #59757-1 50ul #59757-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SULT2A1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	SULT2A1 Antibody detects endogenous levels of total SULT2A1
Immunogen Description	A synthesized peptide derived from human SULT2A1
Other Names	HST; ST2; STD; hSTa; DHEAS; ST2A1; ST2A3; DHEA-ST;
Accession No.	Uniprot:Q06520
Uniprot	Q06520
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:100~1:500 ICC/IF 1:50~1:200 FC 1:80

Images



Western blot analysis of SULT2A1 expression in HepG2 cell lysate.

Product Description

Catalyzes the sulfation of steroids and bile acids in the liver and adrenal glands.

Background

Catalyzes the sulfation of steroids and bile acids in the liver and adrenal glands.

Note: This product is for in vitro research use only