

IRAK2 Rabbit mAb

Catalog No: #59760

Package Size: #59760-1 50ul #59760-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IRAK2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP FC
Species Reactivity	Human
Specificity	IRAK2 Antibody detects endogenous levels of total IRAK2
Immunogen Description	A synthesized peptide derived from human IRAK2
Other Names	Interleukin 1 receptor associated kinase 2; IRAK 2; Irak2;
Accession No.	Uniprot:O43187
Uniprot	O43187
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IP 1:50 FC 1:100

Images



Western blot analysis of IRAK2 expression in HeLa cell lysate.

Product Description

Binds to the IL-1 type I receptor following IL-1 engagement, triggering intracellular signaling cascades leading to transcriptional up-regulation and mRNA stabilization.

Background

Binds to the IL-1 type I receptor following IL-1 engagement, triggering intracellular signaling cascades leading to transcriptional up-regulation and mRNA stabilization.

Note: This product is for in vitro research use only