

## Oct-1 Rabbit mAb

Catalog No: #59761

Package Size: #59761-1 50ul #59761-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Oct-1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	Oct-1 Antibody detects endogenous levels of total 43374
Immunogen Description	A synthesized peptide derived from human 43374
Other Names	NF A1; OCT1; OTF1; PO21; PO2F1; POU2F1;
Accession No.	Uniprot:P14859
Uniprot	P14859
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:1000 IHC 1:50~1:200

## Images



Western blot analysis of 43374 expression in HEK293 cell lysate.

## Product Description

Transcription factor that binds to the octamer motif (5'-ATTTCAT-3') and activates the promoters of the genes for some small nuclear RNAs (snRNA) and of genes such as those for histone H2B and immunoglobulins.

## Background

---

Transcription factor that binds to the octamer motif (5'-ATTTGCAT-3') and activates the promoters of the genes for some small nuclear RNAs (snRNA) and of genes such as those for histone H2B and immunoglobulins.

---

Note: This product is for in vitro research use only