

Cellubrevin Rabbit mAb

Catalog No: #59762

Package Size: #59762-1 50ul #59762-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Cellubrevin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human
Specificity	Cellubrevin Antibody detects endogenous levels of total Cellubrevin
Immunogen Description	A synthesized peptide derived from human Cellubrevin
Other Names	CEB; Cellubrevin; Synaptobrevin 3; VAMP3;
Accession No.	Uniprot:Q15836
Uniprot	Q15836
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:40 FC 1:50

Images



Western blot analysis of Cellubrevin expression in A375 cell lysate.

Product Description

SNARE involved in vesicular transport from the late endosomes to the trans-Golgi network.

Background

SNARE involved in vesicular transport from the late endosomes to the trans-Golgi network.

Note: This product is for in vitro research use only