

FHL2 Rabbit mAb

Catalog No: #59764

Package Size: #59764-1 50ul #59764-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FHL2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	FHL2 Antibody detects endogenous levels of total FHL2
Immunogen Description	A synthesized peptide derived from human FHL2
Other Names	AAG11; DRAL; Fhl2; SLIM3;
Accession No.	Uniprot:Q14192
Uniprot	Q14192
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:1000 ICC/IF 1:50~1:200 IP 1:40

Images



Western blot analysis of FHL2 expression in SW480 cell lysate.

Product Description

May function as a molecular transmitter linking various signaling pathways to transcriptional regulation. Negatively regulates the transcriptional repressor E4F1 and may function in cell growth.

Background

May function as a molecular transmitter linking various signaling pathways to transcriptional regulation. Negatively regulates the transcriptional repressor E4F1 and may function in cell growth.

Note: This product is for in vitro research use only