

p53 DINP1 Rabbit mAb

Catalog No: #59766

Package Size: #59766-1 50ul #59766-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	p53 DINP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	p53 DINP1 Antibody detects endogenous levels of total p53 DINP1
Immunogen Description	A synthesized peptide derived from human p53 DINP1
Other Names	p53DINP1; SIP; Teap; TP53DINP1; TP53INP1; TP53INP1A; TP53INP1B;
Accession No.	Uniprot:Q96A56
Uniprot	Q96A56
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of p53 DINP1 expression in (1) HepG2 cell lysate; (2) RAW 264.7 cell lysate; (3) PC12 cell lysate.

Product Description

In response to double-strand DNA breaks, promotes p53/TP53 phosphorylation on 'Ser-46' and subsequent apoptosis.

Background

In response to double-strand DNA breaks, promotes p53/TP53 phosphorylation on 'Ser-46' and subsequent apoptosis.

Note: This product is for in vitro research use only