

Apolipoprotein A II Rabbit mAb

Catalog No: #59780

Package Size: #59780-1 50ul #59780-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Apolipoprotein A II Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Rat
Specificity	Apolipoprotein A II Antibody detects endogenous levels of total Apolipoprotein A II
Immunogen Description	A synthesized peptide derived from human Apolipoprotein A II
Other Names	APO A2; Apo AII; APOA 2; ApoA II; APOA2; APOAII; Apolipoprotein A II;
Accession No.	Uniprot:P02652
Uniprot	P02652
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of Apolipoprotein A II expression in Human plasma lysate.

Product Description

May stabilize HDL (high density lipoprotein) structure by its association with lipids, and affect the HDL metabolism.

Background

May stabilize HDL (high density lipoprotein) structure by its association with lipids, and affect the HDL metabolism.

Note: This product is for in vitro research use only