

Tetranectin Rabbit mAb

Catalog No: #59785

Package Size: #59785-1 50ul #59785-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Tetranectin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	Tetranectin Antibody detects endogenous levels of total Tetranectin
Immunogen Description	A synthesized peptide derived from human Tetranectin
Other Names	Clec3b; Tetranectin; TN; TNA;
Accession No.	Uniprot:P05452
Uniprot	P05452
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:100~1:500 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of Tetranectin expression in Human plasma lysate.

Product Description

Tetranectin binds to plasminogen and to isolated kringle 4. May be involved in the packaging of molecules destined for exocytosis.

Background

Tetranectin binds to plasminogen and to isolated kringle 4. May be involved in the packaging of molecules destined for exocytosis.

Note: This product is for in vitro research use only