

RanGAP1 Rabbit mAb

Catalog No: #59793

Package Size: #59793-1 50ul #59793-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	RanGAP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	RanGAP1 Antibody detects endogenous levels of total RanGAP1
Immunogen Description	A synthesized peptide derived from human RanGAP1
Other Names	SD; Fug1; RANGAP;
Accession No.	Uniprot:P46060
Uniprot	P46060
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:20 FC 1:50

Images



Western blot analysis of RanGAP1 expression in HeLa cell lysate.

Product Description

GTPase activator for the nuclear Ras-related regulatory protein Ran, converting it to the putatively inactive GDP-bound state.

Background

GTPase activator for the nuclear Ras-related regulatory protein Ran, converting it to the putatively inactive GDP-bound state.

Note: This product is for in vitro research use only